Solutions for safe transfers and mobility



Contents



Safe transfers and handling

- 04 VertifHelp transfer aid
- 04 Etac® Molift Raiser Pro
 - 05 Etac® Molift Raiser Pro accessories
- 06 Etac® Immedia transfer boards.
 - 06 Etac® Immedia 3B Board
 - 07 Etac® Immedia E-board
 - 07 Etac® Immedia Butterfly board with GlideWing
- 08 Etac® Immedia SatinSheets.
 - 08 Etac® Immedia SatinSheet 4Direction bed sheet
 - 09 Etac® Immedia SatinSheet 2Direction base sheet
 - 09 Etac® Immedia SatinSheet 2Direction pillowcase
- 10 Hepro transfer assistance system
- Etac® Molift RgoSling ambulating vest
- 12 Silvalea® slings.
 - 00 Silvalea® Universal deluxe sling
 - 00 Silvalea® Superfine Plus in-situ deluxe sling



Mobility

- 14 Pacer® gutter walking frame.
- 14 Nomad® gutter walking frame
- 15 Pacer® walking frame
- 15 Mobilis® Narrow super stroller-8 walking frame
- 16 Freiheit® Freedom stroller XC walking frame
- 16 Individualise your walking frame
- 17 Cobi Rehab XXL mobility aids
 - 17 Cobi Rehab XXL walking aid
 - 17 Cobi Rehab XXL crutch
 - 18 Cobi Rehab XXL walking frame
- 18 Cobi Rehab XXL Stand tall walking frame
- 19 Cobi Rehab XXL Rollator King
- 0 Super Stroller® Maxi steel walking frame



Cubro

Safe transfers and handling

VertiHelp transfer aid



Code: 67785

Max user weight: 180kg

The VertiHelp transfer aid is perfect for facilitating quick and smooth transfers. It has been designed for people who have difficulty standing up on their own from a seated position and is used to transfer between seated surfaces. The user is lifted to a partial or full standing position while both feet remain on a secure platform.

It has a small turning radius making it ideal for smaller spaces. A non-powered device, simplifying everyday transfers.

Etac® Molift Raiser Pro



Code: 66553 SWL: 170kg

- The ergonomic platform for a safe and active sit-to-stand transfer
- Multiple grips for both the user and carer
- Two soft padded leg supports can be height adiusted
- Low base profile to fit under chairs
- The front of the platform is curved to fit closer to the toilet
- Central brake mechanism is activated by one foot on either side



Padded lea supports tilt independently



Multiple hand grips for both the carer and the

Etac® Molift Raiser Pro ACCESSORIES








Code: S/M 66554

Code: L/XL 66555

Assists the user up to a standing position.







Code: S/M 66559

Code: L/XL 66560

Allows for a smooth transfer when assisting the user up to a standing position.





Code: 66864



Molift Raiser Pro heel strap Molift Raiser Pro block

Code: 66863

Etac[®] Immedia transfer boards

The Immedia range of transfer boards safely bridge the gap when performing seated transfers. The smooth products reduce friction and facilitate weight shifts.

Etac® Immedia 3B Board



Code: 66438

A stable transfer board with a traditional curved shape. The board can be placed either curve forward or curve backward depending on the transfer situation and surrounding obstacles e.g. the wheelchair wheel.

The board's low friction upper surface allows the user to transfer independently or with assistance.







Etac® Immedia E-board



Code: (S) 66439 Code: (M) 66440

The Immedia E-Board is a thin and flexible transfer board that allows easier application and removal. Can be used for transfers from bed to wheelchair or shower chair.

It can also be used for transfers in and out of the car and between wheelchair, chair, toilet, shower chair and commode.

Etac® Immedia Butterfly board with GlideWing



Code: (L) 66443/66444 Code: (XL) 66443/66445

Long, wide and stable transfer board intended for seated transfers. Featuring cut-outs for positioning against the wheel of a wheelchair and non-slip pads on the underside.

The GlideWing slides along the board and provides even lower friction and is recommended when the transfer involves bare skin, e.g. shower chair to toilet.

6 Solutions for safe transfers and mobility 0800 101 211 | communitycare@cubro.co.nz | cubro.co.nz 7



Cubro

Etac[®] Immedia SatinSheets

Designed to make moving an individual in bed easier, whatever their ability. The in-bed system helps reduce strain on the carers back and therefore prevents work related injuries.

As a result, helps save time, may reduce the number of carers required and ensures the individual is comfortable and well supported.

Etac® Immedia SatinSheet 2Direction base sheet

Code: (Single) 66316 Code: (Super single) 66317 Code: (Double) 66318





Elastic corner straps to fit any mattress





View the in-bed system in action

Check out our video demonstrations of how the SatinSheet in-bed positioning system works at cubro.co.nz/immedia

Etac® Immedia SatinSheet 4Direction bed sheet

Code: (Midi) 66322 Code: (Maxi) 66324









(1) Handles for ease of movement

Cubro

Hepro transfer assistance system





Code: 4055 SWL: 150kg

The Hepro transfer assistance system is a sturdy and safe support aid that helps individuals to stand and move around their homes. The telescopic support pole can be used on its own or supplemented with a range of useful accessories to suit any home environment and individual needs.

The Hepro system provides users with greater independence and confidence to perform daily tasks. It fits in every home and is very easy to install. It simply clamps between the floor and ceiling so does not require the use of screws, nails or glue.

ACCESSORIES



Assistance bar Code: 4060



Triangle grip Code: 4061



Hand grip Code: 4062



Ceiling plate

Code: 4070


Etac® Molift RgoSling ambulating vest



SWL: 300kg

Developed for safe and secure rehabilitation such as gait training, standing exercise as well as balance exercise.

- A double buckled wide padded belt that distributes pressure evenly over a larger surface
- The self-tightening suspension loops prevent users from slipping or falling down





Padded groin straps limit the vest from sliding up Codes below include padded groin straps to limit the vest from sliding up. Also available without groin straps.

WAIST SIZE		CODE
XXS	500 - 650mm	66835
XS	600 - 750mm	66836
S	700 - 900mm	66837
М	850 - 1050mm	66838
L	1000 - 1300mm	66839
X	1250 – 1550mm	66840

10 Solutions for safe transfers and mobility communitycare@cubro.co.nz | cubro.co.nz | 11

Silvalea® slings

Silvalea® Universal deluxe sling with head support



SWL: 220kg

Full body sling designed with shaped leg sections providing support to both the hip and thigh.

- Suitable for active users and asymmetrical body shapes
- High density open weave mesh fabric is fast drying, suitable for bathing
- Some access for personal hygiene
- Attachment type: loops



Silvalea® Universal deluxe sling with head support



SWL: 220kg

Provides full body support on transfer and offers hip, thigh and full femur support.

- Detachable tapes for greater comfort in-situ
- High-performance fabric offering superior comfort, suitable for bathing
- Attachment type: loops



WAIST SIZE	CODE
XXXS	6606S
XXS	6605S
XS	6600S
S	6601S
М	6602S
L	6603S
XL	6604S



Cubro

Mobility

Cubro

Pacer® gutter walking frame



With 8" wheels

Code: 66796 SWL: 160ka

The Pacer® gutter walking frame is designed with vertical autter supports resulting in an exceptionally balanced frame. The gutter troughs allow the user to bear weight through their forearms rather than hands which is useful for individuals who have arthritis or an injury to their hand or wrist.

- Height and angle adjustable gutter troughs
- Sturdy steel frame
- Knurled gutter stems prevent slippage and rotation
- "One touch" brake system
- Folds easily for storage and travel
- Comes with basket and tray
- Suited to both indoors and outdoors

Nomad® gutter walking frame



Code: 1885HB SWL: 190kg

The Nomad® gutter walking frame is an exceptionally rigid steel frame with forearm gutter support. Its narrow front assists manoeuvring around furniture and through doorways. A high front bar allows the frame to be wheeled over a toilet.

- Height and angle adjustable gutter troughs
- Suited for indoor use
- Steel frame
- Knurled gutter stems to prevent slippage and rotation
- Folds down for transport



Pulpit top

Code: 1885/1867

Pacer® walking frame



With 8" wheels and handbrakes

Code: 66797 SWL: 160kg

The Pacer® walking frame is exceptionally balanced and is made from high-quality materials with the steel frame built for durability and strength.

- Suitable for indoor and outdoor use, with 8" wheels that make it a lot easier to push and manoeuvre, especially outdoors
- · Sturdy steel frame
- · "One touch" brake system
- Folds easily for storage and travel
- · Height adjustable handles
- Anatomical handgrips
- · Comes with basket and trav

Mobilis® Narrow super stroller-8 walking frame



Code: 66795 SWL: 120kg

This frame is shorter, narrower and lower than most standard walking frames meaning it is perfect for use in space restricted areas and for people of smaller stature. Its larger 200mm (8") wheels make it suitable for indoor and outdoor use.

14 Solutions for safe transfers and mobility 0800 101 211 | communitycare@cubro.co.nz | cubro.co.nz 15

Freiheit® Freedom stroller XC walking frame



Code: 1554 SWL: 150kg

The most robust walking frame in our Freiheit® freedom stroller range, the XC is designed for more adventurous individuals. For those who like to get out and about, the large 310mm wheels and excellent stability make it easy to push the frame over rough ground, making it perfect for outdoor use.

- Large 310mm pneumatic wheels
- Seat with backrest
- Mesh bag holds up to 5kg
- Walking stick holder
- Titling pedal to overcome obstacles
- · Robust tubular steel frame
- Easy to fold and lock for transport
- Reflectors for enhanced safety
- Steering dampers reduce the 'flutter' of the steering wheels on rough terrain

Individualise your walking frame

We can individualise our walking frames to meet all your client's needs. Ask us about:

- · Cubro adjustable pulpit top to suit Cubro Gutter frames
- · Pacer pram handle to suit
- · Pacer walking frames
- Super stroller oxygen bottle holder
- Super stroller folding tray
- Super stroller under seat bag

Cubro

Cobi Rehab XXL mobility aids

Mobility solutions designed to support bariatric users with independence, comfort and rehabilitation.

Cobi Rehab XXL walking aid



Code: 9556 SWL: 270ka The reinforced walking aid is specifically designed for bariatric users. Adjustable in height and the swan neck positions the user's weight directly over the shaft for stability and protects the wrist from strain.

Also available with a 4-point base

Code: 9557

Cobi Rehab XXL crutch



Anatomical

Code: 9556 (Universal grip)

Code: 9556 (Anatomical grip)

SWL: 325kg

The single crutch optimises comfort with forearm support and is adjustable in height from 680 - 1100mm. The universal grip crutch provides flexibility for the user to alternate sides.

← Universal

0800 101 211 | communitycare@cubro.co.nz | cubro.co.nz 17



Cobi Rehab XXL walking frame



Code: 9553 (740mm) Code: 9554 (840mm)

SWL: 270kg

The strona, functional and liahtweight design features a height adjustable frame that can be adjusted from 760mm -1010mm.

- Folds easily for transport and storage
- Padded handles provide the user with added comfort

Cobi Rehab XXL Stand tall walking frame



Code: 9552

The Cobi Rehab XXL Stand Tall walking frame is a multi-functional mobility aid that can be used for walking and sit-to-stand transfers. The Stand Tall gives greater independence as it allows users to engage in their own transfers.

Specifically designed for bariatric users, this mobility aid has a user weight of 325kg. 32mm tubing ensures the frame is durable and long lasting as well as easy to clean.

- · Easily adjustable frame height
- The width and angle of the forearm padding is adiustable
- Lockable brakes
- Durable, easily cleaned frame constructed from 32mm diameter tubing
- · Very high user weight for bariatric users
- 2x 100mm front castors, 2x 125mm rear castors

Cobi Rehab XXL Rollator King



Code: 9555 SWL: 325kg

A reinforced rollator specifically designed for bariatric users. Liahtweight and suitable for indoor and outdoor use, the rollator supports a very high user weight of 325kg.

The Rollator King fits through doorways and folds easily for transport and storage by lifting the seat.

- 4x 200mm wheels
- Highly responsive 'feather touch' hand brake system
- Folds down for transport and storage
- Basket attached
- · Powder-lacquered aluminium frame
- Large 510mm wide supportive seat with back
- Reflectors fitted to handles

Super Stroller® Maxi steel walking frame



Code: 1955

Built to last, the Super Stroller walking frame has a reinforced, heavy-duty, steel constructed frame, with a high user weight. The wide frame, larger footprint and high user weight make it a great option for larger users.

The 200mm wheels have steel vokes and are suitable for both indoor and outdoor use.

- Storage basket included
- · Heavy-duty steel frame with steel yokes on the wheels for added strength
- · Large 200mm front castors and rear wheels,
- · Easy to adjust "one touch" handbrake
- Padded seat
- · Easy to fold for storage

18 Solutions for safe transfers and mobility 0800 101 211 | communitycare@cubro.co.nz | cubro.co.nz 19

